

# BRIGHTEN YOUR OPERATING ROOM WITH RAOUF COMPLETE SYSTEM

## Camera Head

- CMOS sensor with full HD resolution (1920 x 1080 pixels)
- Integrated 16-34mm HD zoom
- 9 Independently programmable control buttons
- Ergonomically designed by physicians and clinicians

Sensor twice as sensitive\*, for great image quality even in low light conditions



## Cam+

### Camera Control Unit:

- Presets programming for LAP, ENT, GYN, URO, Arthroscopy and Flexible Endoscopy
- NFC technology allows individual physician SCard programming
- Unlimited physician profile settings
- Integrated communication with Led and Flow
- Integrated Full HD image and Full HD video recording capability
- All functions are programmable for camera head use
- Touch screen with 16:9 resolution
- Easy to plug monofilament connector



**EXCLUSIVE**

BlueTest function to enhance blue colors for specific usages



**PATENTED**



## Flow

### Insufflator:

- Insufflation can be controlled directly from the camera head
- From 2L/min to 45L/min insufflation range
- Automatic management of insufflation flow to maintain desired pressure
- Proprietary technology to calculate remaining available insufflation time
- Insufflation pressures displayed simultaneously on both Monitor and Touch screen
- Large font, easy to read touch screen with 16:9 resolution
- Internal heating system
- Bariatric specific insufflation parameters
- External desufflation valve to minimize cross contamination and over insufflation
- Great efficiency in low pressure (8-9mmHg)

Remaining insufflation time displayed in minutes



## Led

### LED Light Source:

- 50,000h lifetime
- LED intensity equivalent to a 300W Xenon light source
- Touch screen with 16:9 resolution
- Automatic light cable detection for safety
- Automatic luminosity adjustment
- Light can be controlled directly from the camera head



SONY Medical Monitors

**CONNECTED SOLUTION FOR SAFE SURGERIES**



**COMPLETE SOLUTION**



# DUO

Full HD

The work quality is closely linked to the ability to discern the smallest details. Having control of these details makes the difference!

Good quality of light is essential for any endoscopic act, this is why it is important to allow the practitioner to adapt the parameters according to his requirements and the type of treatment.


Full HD camera

The Full HD resolution provides accurate images to ensure a reliable treatment.

Led Light source



## ubicam



شرکت سونیران  
آدرس: تهران - خیابان سهروردی شمالی - خیابان خرمشهر -  
خیابان عربعلی - کوچه هفتم - پلاک ۲۴ و ۲۵  
تلفن: ۹۱۰۰۹۴۱۴ فکس: ۸۸۷۴۱۴۶۷

DISCOVER THE LATEST  
FULL HD TECHNOLOGY



**SONIRAN**  
MEDICAL GROUP

## HIGHLIGHT YOUR SURGERIES

