

TECHNICAL SPECIFICATIONS

Catheter System	: Over-the-wire (OTW)
Drug Type	: Paclitaxel
Drug Dosage	: 3 µg/mm ²
Usable Catheter Length	: 80 cm & 135 cm
Max. Guidewire Compatibility	: 0.035"
Balloon Diameter (mm)	: 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00 & 10.00 mm
Balloon Length (mm)	: 30, 40, 50, 60, 80 & 120 mm
Min. Sheath Compatibility	: 5 Fr to 8 Fr (Refer product label/IFU for more details)
Balloon RO Markers	: 2 - Platinum/Iridium
Shaft Diameter	: 5 Fr
Nominal Pressure(NP)	: 6-7 atm depending upon diameter and length of balloon (Refer product label/IFU for more details)
Rated Burst Pressure (RBP)	: 11-14 atm depending upon diameter and length of balloon (Refer product label/IFU for more details)

Mozec™ PEB PTA - ORDERING INFORMATION

Usable Catheter length 80 cm

Balloon Diameter (mm)	Balloon Length (mm)					
	30	40	50	60	80	120
3.00	MO35P030030A	MO35P030040A	MO35P030050A	MO35P030060A	MO35P030080A	MO35P030120A
4.00	MO35P040030A	MO35P040040A	MO35P040050A	MO35P040060A	MO35P040080A	MO35P040120A
5.00	MO35P050030A	MO35P050040A	MO35P050050A	MO35P050060A	MO35P050080A	MO35P050120A
6.00	MO35P060030A	MO35P060040A	MO35P060050A	MO35P060060A	MO35P060080A	MO35P060120A
7.00	MO35P070030A	MO35P070040A	MO35P070050A	MO35P070060A	MO35P070080A	MO35P070120A
8.00	MO35P080030A	MO35P080040A	MO35P080050A	MO35P080060A	MO35P080080A	MO35P080120A
9.00	MO35P090030A	MO35P090040A	MO35P090050A	MO35P090060A	MO35P090080A	-
10.00	MO35P100030A	MO35P100040A	MO35P100050A	MO35P100060A	MO35P100080A	-

Usable Catheter length 135 cm

Balloon Diameter (mm)	Balloon Length (mm)					
	30	40	50	60	80	120
3.00	MO35P030030B	MO35P030040B	MO35P030050B	MO35P030060B	MO35P030080B	MO35P030120B
4.00	MO35P040030B	MO35P040040B	MO35P040050B	MO35P040060B	MO35P040080B	MO35P040120B
5.00	MO35P050030B	MO35P050040B	MO35P050050B	MO35P050060B	MO35P050080B	MO35P050120B
6.00	MO35P060030B	MO35P060040B	MO35P060050B	MO35P060060B	MO35P060080B	MO35P060120B
7.00	MO35P070030B	MO35P070040B	MO35P070050B	MO35P070060B	MO35P070080B	MO35P070120B
8.00	MO35P080030B	MO35P080040B	MO35P080050B	MO35P080060B	MO35P080080B	MO35P080120B
9.00	MO35P090030B	MO35P090040B	MO35P090050B	MO35P090060B	MO35P090080B	-
10.00	MO35P100030B	MO35P100040B	MO35P100050B	MO35P100060B	MO35P100080B	-

Sheath Compatibility(OTW 0.035")

5Fr for balloon diameters 3.00, 4.00 and 5.00 mm
 6Fr for balloon diameters 6.00 and 7.00 mm
 7Fr for balloon diameter 8.00 mm
 8Fr for balloon diameters 9.00 and 10.00 mm

Meril

More to Life

Manufacturer:

Meril Life Sciences Pvt. Ltd.
 Survey No. 135/139, Bilakhia House,
 Muktanand Marg, Chala, Vapi - 396 191.
 Gujarat, India.
 T +91 260 240 8000

Subsidiary companies:

Meril Life Sciences Pvt. Ltd.
 301, A-Wing, Business Square,
 Chakala, Andheri Kurla Road,
 Andheri East, Mumbai 400 093
 T +91 22 39350700
 F +91 22 39350777

Meril, Inc.

2436 Emrick Boulevard,
 Bethlehem, PA 18020
 T +610 500 2080
 F +610 317 1672

Meril South America

Meril South America
 Doc Med LTDA
 Alameda dos Maracatins 992,
 Edificio Aquarius, Andar 8,
 Bloco B, Moema - CEP -04089001,
 Sao Paulo, Brazil.
 Tel: +5511 2506-3624

Meril GmbH.

Bornheimer Strasse 135-137,
 D-53119 Bonn,
 Germany.
 T +49 228 7100 4000
 F +49 228 7100 4001

Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş.
 İçerenköy Mah. Çetinkaya Sok.
 Prestij Plaza No:28
 Kat:4 Ataşehir, 34752
 İstanbul / Turkey
 T +90 216 641 44 24
 F +90 216 641 44 25

Meril China Co. Ltd.

2301b,23f, Lixin Plaza,
 no 90, South Hubin Road,
 Xiamen, China
 T 0086-592-5368505
 F 0086-592-5368519

Meril SA Pty. Ltd.

102, 104, S101 and S102,
 Boulevard West Office Park,
 142 Western Service Road,
 Erf 813 Woodmead Extension 17
 Sandton, Johannesburg - 2191
 South Africa
 T +27 11 465-2049
 F +27 86 471 7941

Meril Medical LLC.

Nauchnyi Proezd 19,
 Moscow, Russia - 117 246.
 Office - +7 495 772 7643


EU representative.

Obelis S.A.
 Bd. General Wahis 53,
 1030, Brussels, Belgium.
 T +32 2 732 5954
 F +32 2 732 6003
 E mail@obelis.net

E askinfo@merillife.com
 W www.merillife.com

Meril

More to Life



MOZEC™ PEB PTA

Paclitaxel Eluting PTA Balloon Catheter (OTW 0.035")

Contact your country local Meril sales representative for availability of sizes.
 Mozec™ PEB PTA is a registered trademark of Meril Life Sciences Pvt. Ltd.
 Mozec™ PEB PTA is not approved and not available for sale in USA.

CE
1783

MO35P/IBROCHURE/005/MLS/20211110


CE

MOZEC™ PEB PTA

Paclitaxel Eluting PTA Balloon Catheter (OTW 0.035")

Meril

More to Life



Size Matrix :

Ø : 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00 & 10.00 mm

ℓ : 30, 40, 50, 60, 80 & 120 mm