



# CLISIS *Evolution*

R/F Remote-controlled  
system with digital flat  
panel detector



GMM  
GROUP

# THE EVOLUTION OF A TOP QUALITY PRODUCT

REDESIGNED, REENGINEERED, RESTYLED TO OFFER YOU  
THE BEST EXPERIENCE EVER

## OUTSTANDING PERFORMANCES

CLISIS is the **evolution** of a remote-controlled table, representing the best combination of operational **functionality**, premium **design** and enhanced image **quality** in any radiographic and fluoroscopic procedure

- **Advanced UI**
- **New design**
- Outstanding operational **efficiency**
- Complete **reliability**
- State-of-the-art **digital flat panel detector**
- Powerful **image processing**
- Fully **DICOM** compliant





## TOTAL ACCESSIBILITY

**The rear accessibility** offered by the CLISIS tabletop allows easy patient positioning, granting maximum **comfort** and **efficiency**.



ALL IN ONE:  
FUNCTIONALITY AND DESIGN







- **90/90** remote tilting system
- Wide adjustable **SID**
- **Reduced** tabletop-floor **distance**
- Single-end suspended carbon fiber tabletop
- Enhanced maximum patient **weight**
- Wide displacement of tube/ detector assembly for a **full-length patient coverage**



## WIDE RANGE OF APPLICATIONS

The new CLISIS brings to **higher innovation** and **optimized workflow** both in routine and specialized procedures: continuous and pulsed fluoroscopy, high resolution digital spot exposures, tomography, tomosynthesis, spine and lower limbs automatic image stitching.



OBLIQUE PROJECTION



EXAMINATIONS IN  
TRENDELENBURG POSITION



UROGENTIAL STUDIES



SPECIAL PROJECTIONS  
ON STRETCHER



SPECIAL PROJECTIONS WITH  
WIRELESS DETECTOR



GASTROINTESTINAL  
EXAMINATIONS



WIDE RANGE OF ACCESSORIES

## UTMOST VERSATILITY FOR SPECIALIZED APPLICATIONS

Easy examinations of patients  
with **reduced mobility**



## EXAMINATIONS WITH HEIGHT ADJUSTABLE FOOTREST





## EMERGENCY APPLICATIONS: SAFETY AND AGILITY

A smooth stretchers' transfer ensures **immediate intervention** also in emergency applications.




# INNOVATIVE USER INTERFACE

**New integrated console:**  
complete remote and quick control of the overall system



**Redesigned ergonomics  
and intuitive user interface:**  
maximum efficiency for the operator





## DIAGNOSIS: QUALITY AND ACCURACY

Significantly improved processes, dose reduction and image enhancement is the result of a **renovated workflow** from the image acquisition to its storage and transmission.

Each specific exam with fluoro, rad, tomography and tomosynthesis procedure automatically elaborates processing algorithms, reducing **waiting times** and displaying **optimized images** almost in real time

## MINIMUM PATIENT DOSE

Ideally suited for **pediatric procedures** thanks to the minimum dose delivered and the possibility to work in **direct contact** with the detector






CLUSIS EVOLUTION EN - 5565 - 0121 Data is subject to change without prior notice

## GET IN TOUCH

GENERAL MEDICAL MERATE SPA  
Via Partigiani, 25  
24068 SERIATE (BG) - Italy

Phone +39 035 45 25 311  
Fax +39 035 29 77 87  
info@gmmspa.com



[www.gmmspa.com](http://www.gmmspa.com)